

Raman spectroscopy

Leena Pitkänen
28.4.2025

1

Raman spectroscopy?

Spectroscopic technique to determine vibrational modes of molecules



Structural information


- Chemical structure, functional groups
- Crystallinity
- Molecular interactions



C. V. Raman

2


Electromagnetic radiation



- Electromagnetic radiation in the near visible light changes the status of molecular vibrations
 - Vibrational spectra of molecules


3

Behavior of light



4

Light scattering



Light scattering can be either elastic or inelastic:

1. Elastic scattering: scattered light has the **same** frequency than incident radiation (Rayleigh).
2. Inelastic scattering: scattered light has **higher** frequency than the incident radiation (anti-Stokes).
3. Inelastic scattering: scattered light has **lower** frequency than the incident radiation (Stokes → Raman).

5

Rayleigh scattering

- Most of the scattered light is at the same wavelength as the incident light

➡ Rayleigh scattering


Useful when measuring polymer
molar masses in solution!

$$LS \sim M \cdot c \cdot (\partial n / \partial c)^2$$

6

Inelastic scattering

The photon frequency decreases
→ the final energy level of molecular vibration is **higher** than the initial energy



The photon frequency increases
→ the final energy level of molecular vibration is **lower** than the initial energy

Stronger intensity

7

Raman (Stokes) scattering

- **Small amount of light (typically 0.0000001%) is scattered at different wavelength than the incident light**
 - This scattering depends on the chemical structure of the analyte
- Raman spectrum commonly shows a number of peaks showing the intensity and wavelength position of the Raman scattered light
- Each peak corresponds to a specific molecular bond vibration

C – C

C – H


C = C

N – O

8

Molecular vibrations

Normal modes of vibrations in a CO₂ molecule




Can be detected with Raman

9

Raman active molecules

- In Raman active molecules the vibration mode must cause polarizability changes
 - Change in size, shape or orientation of the electron cloud that surrounds the molecule

The change in ellipsoid electron cloud should be changing during vibration (not return back)




10

Wave number ($\bar{\nu}$)

- Reciprocal of the wavelength in the units of cm^{-1}

$$\bar{\nu} = \frac{1}{\lambda}$$

Absolute wavenumber



Polystyrene raman spectrum | Raman for life
(ramanlife.com)

11

Relative wavenumber

$$\text{Raman wavenumber } [\text{cm}^{-1}] = \frac{10^7}{\lambda_{exc} [\text{in nm}]} - \frac{10^7}{\lambda_{sc} [\text{in nm}]}$$

Example: The 532 nm laser is used and the analyte scatters light at 554 nm. What is the Raman wavenumber?

746 cm^{-1}

12

12

Common Raman vibrations

Functional Group/ Vibration	Region	Raman	InfraRed
Lattice vibrations in crystals, LA modes	10 - 200 cm^{-1}	strong	strong
$\delta(\text{CC})$ aliphatic chains	250 - 400 cm^{-1}	strong	weak
$\nu(\text{Se-Se})$	290 - 330 cm^{-1}	strong	weak
$\nu(\text{S-S})$	430 - 550 cm^{-1}	strong	weak
$\nu(\text{Si-O-Si})$	450 - 550 cm^{-1}	strong	weak
$\nu(\text{Xmetal-O})$	150-450 cm^{-1}	strong	med-weak
$\nu(\text{C-I})$	480 - 660 cm^{-1}	strong	strong
$\nu(\text{C-Br})$	500 - 700 cm^{-1}	strong	strong
$\nu(\text{C-Cl})$	550 - 800 cm^{-1}	strong	strong
$\nu(\text{C-S})$ aliphatic	630 - 790 cm^{-1}	strong	medium
$\nu(\text{C-S})$ aromatic	1080 - 1100 cm^{-1}	strong	medium
$\nu(\text{O-O})$	845 - 900 cm^{-1}	strong	weak
$\nu(\text{C-O-C})$	800 - 970 cm^{-1}	medium	weak
$\nu(\text{C-O-C})$ asym	1060 - 1150 cm^{-1}	weak	strong

13



Common Raman vibrations

$\nu(\text{CC})$ alicyclic, aliphatic chain vibrations	600 - 1300 cm^{-1}	medium	Medium
$\nu(\text{C=S})$	1000 - 1250 cm^{-1}	strong	weak
$\nu(\text{CC})$ aromatic ring chain vibrations	*1580, 1600 cm^{-1}	strong	medium
	*1450, 1500 cm^{-1}	medium	medium
	*1000 cm^{-1}	strong/medium	weak
$\delta(\text{CH}_3)$	1380 cm^{-1}	medium	strong
$\delta(\text{CH}_2)$			
$\delta(\text{CH}_3)$ asym	1400 - 1470 cm^{-1}	medium	medium
$\nu(\text{C-(NO}_2\text{)})$	1340 - 1380 cm^{-1}	strong	medium
$\nu(\text{C-(NO}_2\text{)})$ asym	1530 - 1590 cm^{-1}	medium	strong
$\nu(\text{N=N})$ aromatic	1410 - 1440 cm^{-1}	medium	-
$\nu(\text{N=N})$ aliphatic	1550 - 1580 cm^{-1}	medium	-
$\delta(\text{H}_2\text{O})$	~1640 cm^{-1}	weak broad	strong
$\nu(\text{C=N})$	1610 - 1680 cm^{-1}	strong	medium
$\nu(\text{C=C})$	1500 - 1900 cm^{-1}	strong	weak
$\nu(\text{C=O})$	1680 - 1820 cm^{-1}	medium	strong
$\nu(\text{C}\equiv\text{C})$	2100 - 2250 cm^{-1}	strong	weak
$\nu(\text{C}\equiv\text{N})$	2220 - 2255 cm^{-1}	medium	strong
$\nu(\text{-S-H})$	2550 - 2600 cm^{-1}	strong	weak
$\nu(\text{C-H})$	2800 - 3000 cm^{-1}	strong	strong
$\nu(\text{=C-H})$	3000 - 3100 cm^{-1}	strong	medium
$\nu(\text{=C-H})$	3300 cm^{-1}	weak	strong
$\nu(\text{N-H})$	3300 - 3500 cm^{-1}	medium	medium
$\nu(\text{O-H})$	3100 - 3650 cm^{-1}	weak	strong

14

Raman instruments


- Raman microscopy
- Raman with probe options
- Portable Raman devices



15

Diffraction gratings


- **Disperses Raman-scattered light according to its wavenumbers (before detector)**
 - Diffracts the different wavelengths in discrete angles
- **Gratings affect the spectral resolution (*i.e.* how small wavenumber differences can be detected)**
- **Grating surface contains parallel lines**
 - Each grating has line density l/mm (lines/mm)
 - The higher the line density, the better the spectral resolution
 - With a lower l/mm grating more spectral range will be acquired



16

CCD detector



- **CCD (charged-coupled device) detectors most commonly used**
 - Converts photon signals to electric signals
 - Can detect and record light intensity of discrete wavelengths separated by the diffraction grating



17

Laser wavelengths in Raman

- Raman scattering intensity is proportional to $1/\lambda^4$
 → Raman wavelength ↑
 Raman scattering intensity ↓
- Common wavelengths: 244 nm, 488 nm, 532 nm, 785 nm, 1064 nm




Fewer molecules
absorb in the
near-infrared
region

18

What kind of samples can be analyzed with Raman?

- **Basically anything!**



Pixabay

19

Challenging in respect of Raman...

- **Metallic materials**
 - High reflection of electromagnetic waves
- **Fluorescent samples**
 - *E.g.* lignin
 - Try different wavelengths, time-gated Raman
- **Very dark samples which absorb light**
 - Weak scattering
- **Very polar molecules (H₂O)**
 - Benefit: analytes dispersed in water can be analyzed

20

Comparison of Raman and IR

- Often complementary techniques, both are fast
- Raman works for aqueous samples (IR suffers water absorption effects)
- IR might work better for samples with fluorescence or if the concentration of the analyte is low
- Raman works better in the lower frequencies → crystal lattice vibrations for distinguishing different crystalline forms

21

Sample preparation

- Basically, no sample preparation is need
- The smoother the sample surface, the better spectrum
- Sometimes powder samples or very “fluffy” samples might be challenging

→ compressing the samples



22

Raman is non-destructive method

- In principle, Raman is non-destructive
- Control the laser power in order to avoid any structural changes caused by the heat
 - Samples might burn!
 - Rotation of the sample



23

Is Raman surface active method?

- With confocal Raman it is possible to do depth analysis
 - Depth (Z) resolution
 - Possible to analyze different layers of the sample





FIG. 5. Raman spectra for the PS/PMMA laminate using a dry objective, plotted as a function of apparent depth below the surface.

24

Examples Raman spectra for diamond and polystyrene

Regular array of identical atoms → one peak (one molecular environment)


- Heavy atoms and weak bonds have low Raman shifts
- Light atoms and strong bonds have high Raman shifts



25

Examples Cellulose and modified celluloses


Whole spectrum



26

Examples Cellulose and modified celluloses


Fingerprint region



27


Examples Cellulose and modified celluloses

CH-OH region




28

Examples Carbon materials




29

Examples Carbon materials



30

Examples Carbon materials –carbon nanotubes




A! Aalto University
School of Chemical
Engineering

RBM, radial breathing mode bands → tube expansion and contraction

22.4.2025
31

31


Examples Carbon materials –carbon nanotubes



22.4.2025
32

32


Raman for quality control



33

Raman imaging

- Combines Raman and microscopy to produce “chemical image”
- Raman spectrum from different locations of the sample is collected to build a map in which every pixel consists of an individual spectrum
- The amount of spectral data defines the resolution (commonly less than 1 μm resolution can be achieved)




34

Raman imaging



https://youtu.be/1AhK_3RGAq4



<https://blue-scientific.com/large-area-raman-mapping-fast-high-resolution/>

35


Raman spectroscopy at CHEM

- Renishaw inVia confocal Raman
- Timegate Raman
- Renishaw UV-Raman

36

Renishaw inVia confocal Raman microscope



- Two lasers: 532 nm (green) and 785 nm (NIR)
- Single spectrum acquisition, depth profiles and area mapping with lateral resolution of $< 1\mu\text{m}$
- Spectral resolution of 1 cm^{-1} with high resolution gratings
- Objectives: 20x, 100x (air), 64x (water-immersion)
- Gratings: 2400/1800 l/mm ("high-resolution" = $1\text{ cm}^{-1}\text{ res.}$) and 830 (possibility to measure broader wavenumber range)



37

Timegate Raman


- Allows time resolved measurements
- 532 nm (green) pulsed laser
- Possibility to use probes (BWTek standard probe and touch Raman immersion probe or microscope)



38

Supression of fluorencence

- **Picosecond range laser excitation source and time-gated single photon counting array detector**
 - *Capturing of instantaneous Raman scattering signal while rejecting the longer average delay fluorescence interference*




Renishaw UV-Raman

- **Wavelengths of 244 nm or 257 nm**
- **High excitation with low wavelength laser (efficiency of Raman scattering $1/\lambda^4$)**
- **Less issues with fluorescence (fluorescence typically occurs at wavelengths longer than 300 nm)**




Quantification of lignin and hexenuronic acid using UV-Raman



41

Surface-enhanced Raman spectroscopy (SERS)



- Valence electrons of metal → localized plasmons which enhance incident light as well as scattered signal

42

Raman spectroscopy

Laboratory demos

PUU1 building (Vuorimiehentie 1)

Abio hall, lab 260, 2nd floor

Meet at the lobby!

